



# Award Winning Performance

## DS-11 Series

### Spectrophotometer/Fluorometer

- Full Spectrum UV-Vis
- Integrated Fluorescence
- Powerful Sample QC
- Flexible Export, Wi-Fi, Email, Ethernet, USB
- No Calibration required with SmartPath Technology®



### FL Kits for DNA/RNA Quantification

- Widest Dynamic Range
  - dsDNA 0.5pg/μL -4,000ng/μL
  - RNA 250pg/μL -1,000ng/μL
- dsDNA 20X Greater Sensitivity than competitor Q
- One RNA kit replaces full competitor Q Range

	DS-11 FX+	DS-11 FX	DS-11+	DS-11
Microvolume				
Cuvette				
Fluorescence				

**Request DEMO**

Optional USB fluorescence module  
Available in **RED**, **BLUE**, **WHITE** & **GREY**