



FroggaBio

Hepta-Foot CRYOTUBES

1.8mL Capped Sterile Cryotubes

- Temperature stability of up to -196°C in vapour-phase liquid nitrogen
- 100 Tubes/Bag x 10 bags
- Externally threaded to avoid contamination
- V-bottom for maximum sample recovery
- Autoclavable 121°C for 15 mins
- Hepta-foot structure for better stability
- Certified DNase/RNase/Endotoxin-free

cat# **Only**
CRYO-1.8 **\$199**



FREE SHIPPING
ON ORDERS OVER \$250

Prices are valid until December 31, 2023. Price changes and/or early termination of this promotion may occur at any time. Terms and conditions apply. Free shipping only offered for orders over \$250 CAD. All orders are processed as quickly as possible; if we are experiencing a high volume of orders, shipments may be slightly delayed.



FroggBio

Biocision



CoolCell® LX



Cell Freezing Containers

- Controlled-rate -1°C/minute freezing
- Holds 12 x 1-2mL cryotubes
- No alcohol or fluids required
- Easy & safe to open when frozen
- Proven for stem cells, primary cells, PBMC, cell lines & more!

ALSO AVAILABLE:

NutriFreez® D10
Cryopreservation Media

NutriFreez® D5 Salt-Based
Cryopreservation Media

- Serum & animal component-free
- Referenced over 300 times
- Manufactured under cGMP
- Optimized for maximum cell recovery & protection
- Suitable for various applications, from basic research, through cell banking, to cell therapy



cat#
BCS-405G,
BCS-405O,
BCS-405,
BCS-405PK

Only
\$350