





**Buy 4  
Get 1  
FREE**

# REDSAFE™

- Ethidium bromide (EtBr) alternative
- 20,000X concentrated, just add 5 $\mu$ L to a 100mL agarose gel
- Excitation: 309nm & 419nm, Emission: 537nm
- Detect both dsDNA & RNA



**FREE SHIPPING**

ON ORDERS OVER \$250

Prices are valid until June 30th, 2023. Price changes and/or early termination of this promotion may occur at any time. Terms and conditions apply. Free shipping only offered for orders over \$250 CAD. All orders are processed as quickly as possible; if we are experiencing a high volume of orders, shipments may be slightly delayed.