



1% TAE Agarose Gel

Store at **25°C**
Room Temperature

DNA Mass
(ng/5µL)

Base Pairs

28

10,000

28

8,000

28

6,000

28

5,000

18

4,000

92

3,000

34

2,500

34

2,000

20

1,500

92

1,000

23

750

30

500

45

250

DNA Ladder

11Kb RTU

500µL

From only
\$19.95
per ladder

1 Ladder \$59
cat# DM010-R500

5 Ladders \$179
cat# DM010-R500K \$35.80 ea

10 Ladders \$279
cat# DM010-R500M \$27.90 ea

20 Ladders \$399
cat# DM010-R500G \$19.95 ea

FREE SHIPPING
ON ORDERS OVER \$500

Prices are valid until June 30th, 2023. Price changes and/or early termination of promotions may occur at any time. Free shipping only offered for orders over \$500 USD. Terms and conditions apply.



FroggqBio



ALSO AVAILABLE

Froggarose™
cat# A87-500G

Molecular-grade agarose suitable for routine protein & nucleic acid electrophoresis.

Only
\$350