



CENTRIFUGE TUBES

\$89 ~~\$95~~
/EA

Per 500/Case (Bulk)

When You Order 4 Cases

MIX & MATCH

Cat#: TB50-500, TB15-500

\$119 /EA

Per 500/Case

Racked

Cat#: TR50-500, TR15-500



- ⊙ Polypropylene conical centrifuge tubes
- ⊙ Sterile, DNase/RNase-free, non-pyrogenic
- ⊙ Each case contains 500 tubes (25/pack, TR15-500 contains 50/pack)
- ⊙ Leak-proof seal and accurate graduates
- ⊙ Fast spins and extreme temp: -80°C to 121°C

STERILEEO ②

Cat#	50mL	15mL	Rack	Bulk	Max RCF
TR50-500	✓		✓		9,500 x g
TB50-500	✓			✓	9,500 x g
TR15-500		✓	✓		8,500 x g
TB15-500		✓		✓	8,500 x g