



IM-300D

Digital Microscope for Cell Culture

High-resolution images, 4K camera & 15.6" touch screen.



IM-300D Key Features

Observation modes

- Brightfield and phase contrast.

3 Plan-Achromat phase contrast, LWD objectives

- 10X, 20X, 40X (Quintuple nosepiece).

LED Illumination

- OPTIKA's exclusive X-LED8 light source doubles your light intensity for unmatched performance.

Digital

- 6MP CMOS Camera with USB 3.0; Intel micro PC Windows 11 OS; 15.6" touch screen.

Software

- ProView, for live imaging, acquisitions & analysis.

Accessories

- Customize your system with Optika's range of objectives, sample holders & a mechanical stage.

Request a Quote

Cat#:IM-300D