




# MICROCENTRIFUGE TUBES

**\$29** FOR  
**1,000** TUBES  
WHEN YOU ORDER **10** CASES



- 1.7mL Microcentrifuge tubes
- Easy seal design
- Transparent surface
- Certified RNase, DNase, DNA & Pyrogen free
- Can be centrifuged up to 20,000 RCF
- Temperatures from -80°C to +121°C

Cat#: LMCT1.7B \$39/1000 tubes,  
LMCT1.7B10 \$290/10 cases

Cat#	Description	QTY	Price
LMCT0.5B	0.5 mL Microcentrifuge Tubes	1,000 Tubes	\$29
LMCT1.7B	1.7 mL Microcentrifuge Tubes	1,000 Tubes	\$39
LMCT2.0B	2.0 mL Microcentrifuge Tubes	500 Tubes	\$29
LMCT5.0B	5.0 mL Microcentrifuge Tubes	200 Tubes	\$69