



QuantGene9600 Real-Time PCR System

~~\$25,000~~
\$19,990

- Large Touchscreen Interface
- 96 X 0.2mL Block
- 5-Channel Detection
- High-Temperature Uniformity
- Fast Ramping Speed
- Mobile App Compatible
- Perfect for Qualitative or Quantitative Detection, SNP Analysis, Genotyping, HRM & More



CAT#: FQD-96C

Channel	1	2	3	4	5
Dye/Probe	FAM, SYBR Green 1	HEX/VIC/TET/ JOE/Cy3/ TAMRA	ROX, TEXAS-RED	Cy5, Quasar-670	Cy5-5, Quasar-705