



Syringe Filters

\$69 BUY 4 GET 1
FREE

ANY ASSORTMENT \$55.20 EACH



- Polyethersulfone (PES) syringe filters
- Polypropylene (PP) casing
- Provides the fastest flow rates, low protein binding and low extractables
- For use with tissue culture media, buffers, and biological sample preparations
- Gamma sterilized, DNase, RNase, pyrogen-free
- Individually packaged filters
- Designed with female luer-lock inlet and male luer slip outlet

Cat#	Description	Housing	Price
SF0.22PES	0.22µm PES Syringe Filter, 25mm - 50 per case	PP	\$69
SF0.45PES	0.45µm PES Syringe Filter, 25mm - 50 per case	PP	\$69