

CENTRIFUGE TUBES



- Polypropylene conical centrifuge tubes
- Each case contains 500 tubes (25/pack, TR15-500 contains 50/pack)
- Sterile, DNase/RNase-free, non-pyrogenic
- Leak-proof seal and accurate graduations
- Fast spins and extreme temp: -80°C to 120°C

STERILE EO 

Cat #	50ml	15ml	Rack	Bulk	Max RCF
TR50-500	✓		✓		9,500 x g
TB50-500	✓			✓	9,500 x g
TR15-500		✓	✓		8,500 x g
TB15-500		✓		✓	8,500 x g



FroggaBio
scientific solutions

\$69

per 500/case (Bulk)
when you order 4 cases



Mix and Match

Cat# TB50-500, TB15-500

\$79

per 500/case



Rack / Bulk



FREE SHIPPING
with orders over \$250

Cat#TR50-500, TB50-500,
TR15-500, TB15-500

FREE SHIPPING 
orders over \$250

Prices are valid until June 30th, 2018. Price changes and or early termination of this promo may occur at any time.

230 Canarctic Drive, Toronto, Ontario M3J 2X8, Canada
tel: 416 736 8325 toll free: 877 318 7277 fax: 416 645 3399
e-mail: info@froggabio.com web: www.froggabio.com