



FrogaBio

scientific solutions

FREE SHIPPING orders over \$200

# Syringe Filters

**\$75** Buy 21 **FREE** Get 1

Any Assortment (\$50 each)



- Polyethersulfone (PES) syringe filters
- Polypropylene (PP) casing
- Provide the fastest flow rates, low protein binding and low extractables
- For use with tissue culture media, buffers, and biological sample preparation
- Gamma sterilized, DNase, RNase, pyrogen free
- Individually packaged filters
- Designed with female luer-lock inlet and male luer slip outlet

Cat#	Description	Housing	Price
SF0.22PES	0.22 $\mu$ m PES syringe filter, 25mm –50per case	PP	\$75
SF0.45PES	0.45 $\mu$ m PES syringe filter, 25mm –50per case	PP	\$75