



FroggaBio
scientific solutions

FREE
SHIPPING orders over \$250


Membrane for Western Blot

\$199

Limited Introductory Offer

(35% off, Regular price: \$299)

- Superior strength - will not curl, bend or crack
- Strip and re-probe
- High sensitivity, low backgrounds
- BSA binding capacity up to 100 mg/cm² (NC), 125 mg/cm² (PVDF)
- For use in: WB, Northern, Southern, colony/plaque lift etc



Cat#	Description	Size	Price
TM300-0.22PVDF	PVDF 0.22um	30cm x 3m roll	\$199
TM300-0.45PVDF	PVDF 0.45um	30cm x 3m roll	\$199
TM300-0.22NC	Supported NC 0.22um	30cm x 3m roll	\$199
TM300-0.45NC	Supported NC 0.45um	30cm x 3m roll	\$199